


Rock-120









Multifunctional precision sampling probe

- External & internal mirror polish, external & internal probe washing
- Dedicated sampling probe equipped with sensitive liquid sensor, timely feedback of reagent & sample residuals
- Collision protection, automatic probe depth adjusting



Excellent mixer design

- Teflon coating mixer. No water dropping(reducing carry-over)
- Excellent mixing effect with standard mixing prosedure





Reagent tray & Sample tray

- 80 reagent positions, 40 sample positions
- Non-stop cooling system with peltier pad inside, 24 hours 2°C-14°C
- Bar code reader(optional)



Reaction tray


- 44 reaction cuvettes
- Reaction volume as low as 150µl
- Stable & accurate temperature (37±0.1°C) for reaction




Advanced operating system

- User friendly software interface simple, easy-to-operation
- Sample automatic dilution(decrease, increase, normal)
- Reaction result automatic checking and warning, calibaration result automatic checking, warning alarm log
- Linear & nonlinear calibration, Multi-points calibration up to 8 points
- Quality control apply to westgard, L-T plot, cumulative









 $\textbf{Head Office:} \ \textbf{Building No.6105/6-Quds St, Mokattam, Cairo, EGYPT.}$

Tel: +202 2504 4134 / +201 0001 18003 Fax: +202 2504 4135

Alex Office: Building No.4 Street No.7 From El Naser St. Ezbt Nasr Aldeen , Alexandria, EGYPT.

Tel: +203 4260557 / +201159887706

Website: www.mira-lab.com E-Mail: Sales@mira-lab.com

