



Smart Casa




- Auto Motility
- Auto Morphology
- DNA Fragmentation

Smart Casa Interface:




Sperm Motility :



Main Modules :



Semen Report



Smart Casa Device Components

- 1- MicroScope With Lenses 10X-20X-40X-100X
- 2- Camera HD High Quality with Clear Pictures
- 3- LabTop or Computer All in one
- 4- OS Win (Win 7 or 8 or 10)
- 5- Laser Printer With Colors
- 6- Smart Casa SoftWare Own to Smart Medical

Technical Parameters

- * Maximum of Tested Sperms :1000
- * The Range of Test Speed :0-180 um/s
- * The Frame Number of Picture: 0 -75
- * Resolution of Particle Diameter : 0-150 um/s
- * The Analysis time 1-5 seconds or Longer

Smart Medical For medical and Laboratories Instruments

7 TH Floor - Star Center - Banha - Egypt

smart_medical@yahoo.com

002-01117255260

Qr

